

# SLC7A4 Antibody



PACO05086

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

SLC7A4

**Uniprot**

O43246

**Synonyms:**

SLC7A4; Cationic amino acid transporter 4; CAT-4; CAT4; Solute carrier family 7 member 4

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human CAT-4.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A