

OR2B8P Antibody



PACO05091

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2B8P

Uniprot

P59922

Synonyms:

OR2B8P; OR2B8; Putative olfactory receptor 2B8; Hs6M1-29P

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 2B8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A