

PHACTR4 Antibody



PACO05093

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000,
IHC:1:200-1:1000

Protein Background:

Gene ID:

PHACTR4

Uniprot

Q8IZ21

Synonyms:

PHACTR4; Phosphatase and actin regulator 4

Immunogen:

Synthesized peptide derived from the Internal region of human PHACTR4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A