

Phospho-SNCA (Y133) Antibody



PACO05100

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

SNCA

Uniprot

P37840

Synonyms:

SNCA; NACP; PARK1; alpha-synuclein; Non-A beta component of AD amyloid; Non-A4 component of amyloid precursor; NACP

Immunogen:

Synthesized peptide derived from human Synuclein- α around the phosphorylation site of Y133.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A