

# Phospho-BRCA1 (S1524) Antibody



PACO05116

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:10000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

BRCA1

**Uniprot**

P38398

**Synonyms:**

BRCA1; RNF53; Breast cancer type 1 susceptibility protein; RING finger protein 53

**Immunogen:**

Synthesized peptide derived from human BRCA1 around the phosphorylation site of S1524.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A