

STAT6 Antibody



PACO05143

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

STAT6

Uniprot

P42226

Synonyms:

STAT6; Signal transducer and activator of transcription 6; IL-4 Stat

Immunogen:

Synthesized peptide derived from human Stat6 around the non-phosphorylation site of Y641.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A