

TAT Antibody



PACO05173

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

TAT

Uniprot

P17735

Synonyms:

TAT; Tyrosine aminotransferase; TAT; L-tyrosine:2-oxoglutarate aminotransferase

Immunogen:

Synthesized peptide derived from the Internal region of human TAT.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A