

Phospho-TK1 (S13) Antibody



PACO05212

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:**Gene ID:**

TK1

Uniprot

P04183

Synonyms:

TK1; Thymidine kinase; cytosolic

Immunogen:

Synthesized peptide derived from human Thymidine Kinase around the phosphorylation site of S13.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A