

TK1 Antibody



PACO05213

Product Information

Size:

50ug

Reactivity:

Human, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000

Protein Background:

Gene ID:

TK1

Uniprot

P04183

Synonyms:

TK1; Thymidine kinase; cytosolic

Immunogen:

Synthesized peptide derived from the C-terminal region of human Thymidine Kinase.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A