

TK1 Antibody



PACO05214

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

TK1

Uniprot

P04183

Synonyms:

TK1; Thymidine kinase; cytosolic

Immunogen:

Synthesized peptide derived from human Thymidine Kinase around the non-phosphorylation site of S13.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A