

# TNNI3 Antibody



PACO05223

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000,  
IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

TNNI3

**Uniprot**

P19429

**Synonyms:**

TNNI3; TNNC1; Troponin I; cardiac muscle; Cardiac troponin I

**Immunogen:**

Synthesized peptide derived from human Troponin I-C around the non-phosphorylation site of T142.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A