

TOP2B Antibody



PACO05231

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:**Gene ID:**

TOP2B

Uniprot

Q02880

Synonyms:

TOP2B; DNA topoisomerase 2-beta; DNA topoisomerase II; beta isozyme

Immunogen:

Synthesized peptide derived from the N-terminal region of human Topo II β .

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A