

C3AR1 Antibody



PACO05252

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

C3AR1

Uniprot

Q16581

Synonyms:

C3AR1; AZ3B; C3R1; HNFAG09; C3a anaphylatoxin chemotactic receptor; C3AR; C3a-R

Immunogen:

Synthesized peptide derived from the Internal region of human C3aR.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A