

TRHR Antibody



PACO05253

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat, Monkey

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:**Gene ID:**

TRHR

Uniprot

P34981

Synonyms:

TRHR; Thyrotropin-releasing hormone receptor; TRH-R; Thyroliberin receptor

Immunogen:

Synthesized peptide derived from the Internal region of human TRH-R1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A