

RGPD8/RGPD4/RGPD3 Antibody



PACO05273

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000

Protein Background:**Gene ID:**

RGPD8/RGPD4/RGPD3

Uniprot

O14715/Q53T03/Q7Z3J3/A6NKT7/P0C839

Synonyms:

RANBP2-like and GRIP domain-containing protein 8

Immunogen:

Synthesized peptide derived from the C-terminal region of human RGPD2/3/4/6/8.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A