

# UBA52 Antibody



PACO05274

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

UBA52

**Uniprot**

P62987

**Synonyms:**

UBA52; UBCEP2; Ubiquitin-60S ribosomal protein L40; CEP52; Ubiquitin A-52 residue ribosomal protein fusion product 1

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L40.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A