

COL14A1 Antibody



PACO05282

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:

Gene ID:

COL14A1

Uniprot

Q05707

Synonyms:

COL14A1; UND; Collagen alpha-1(XIV) chain; Undulin

Immunogen:

Synthesized peptide derived from the N-terminal region of human COL14A1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A