

Phospho-VASP (T278) Antibody



PACO05286

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

VASP

Uniprot

P50552

Synonyms:

VASP; Vasodilator-stimulated phosphoprotein; VASP

Immunogen:

Synthesized peptide derived from human VASP around the phosphorylation site of T278.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A