

VAV1 Antibody



PACO05290

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300

Protein Background:

Gene ID:

VAV1

Uniprot

P15498

Synonyms:

VAV1; VAV; Proto-oncogene vav

Immunogen:

Synthesized peptide derived from human Vav around the non-phosphorylation site of Y174.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A