

VTN Antibody



PACO05306

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

VTN

Uniprot

P04004

Synonyms:

VTN; Vitronectin; VN; S-protein; Serum-spreading factor; V75

Immunogen:

Synthesized peptide derived from the Internal region of human Vitronectin.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A