

XRCC3 Antibody



PACO05312

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

XRCC3

Uniprot

O43542

Synonyms:

XRCC3; DNA repair protein XRCC3; X-ray repair cross-complementing protein 3

Immunogen:

Synthesized peptide derived from the Internal region of human XRCC3.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A