

ZNF76 Antibody



PACO05334

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:

Gene ID:

ZNF76

Uniprot

P36508

Synonyms:

ZNF76; D6S229E; ZNF523; Zinc finger protein 76; Zinc finger protein 523

Immunogen:

Synthesized peptide derived from the Internal region of human ZNF76.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A