

SLC30A2 Antibody



PACO05339

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:5000, IHC:1:100-1:300

Protein Background:

Gene ID:

SLC30A2

Uniprot

Q9BRI3

Synonyms:

SLC30A2; ZNT2; Zinc transporter 2; ZnT-2; Solute carrier family 30 member 2

Immunogen:

Synthesized peptide derived from the Internal region of human ZnT-2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A