

SLC30A4 Antibody



PACO05340

Product Information

Size: 50ug	Protein Background:
Reactivity: Human, Mouse, Rat	Gene ID: SLC30A4
Source: Rabbit	Uniprot O14863
Isotype: IgG	Synonyms: SLC30A4; ZNT4; Zinc transporter 4; ZnT-4; Solute carrier family 30 member 4
Applications: ELISA, WB, IHC	Immunogen: Synthesized peptide derived from the Internal region of human ZnT-4.
Recommended dilutions: ELISA:1:5000, WB:1:500-1:2000, IHC:1:100-1:300	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A