

# IRX3 Antibody



PACO05355

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

IRX3

**Uniprot**

P78415

**Synonyms:**

IRX3; IRXB1; Iroquois-class homeodomain protein IRX-3; Homeodomain protein IRXB1; Iroquois homeobox protein 3

**Immunogen:**

Synthesized peptide derived from the Internal region of human IRX3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A