

# IRX1 Antibody



PACO05356

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

ELISA:1:40000, WB:1:500-1:2000,  
IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

IRX1

**Uniprot**

P78414

**Synonyms:**

IRX1; IRXA1; Iroquois-class homeodomain protein IRX-1; Homeodomain protein IRXA1;  
Iroquois homeobox protein 1

**Immunogen:**

Synthesized peptide derived from the Internal region of human IRX1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A