

# OR52N1 Antibody



PACO05361

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR52N1

**Uniprot**

Q8NH53

**Synonyms:**

OR52N1; Olfactory receptor 52N1; Olfactory receptor OR11-61

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 52N1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A