

OR6J1 Antibody



PACO05362

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR6J1

Uniprot

Q8NGC5

Synonyms:

OR6J1; OR6J2; Olfactory receptor 6J1; Olfactory receptor 6J2

Immunogen:

Synthesized peptide derived from the Internal region of human Olfactory receptor 6J1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A