

KLHL36 Antibody



PACO05377

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

KLHL36

Uniprot

Q8N4N3

Synonyms:

KLHL36; C16orf44; Kelch-like protein 36

Immunogen:

Synthesized peptide derived from the C-terminal region of human C16orf44.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A