

EPHX3 Antibody



PACO05380

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:40000, WB:1:500-1:2000,
IHC:1:100-1:300

Protein Background:

Gene ID:

EPHX3

Uniprot

Q9H6B9

Synonyms:

EPHX3; ABHD9; Epoxide hydrolase 3; Abhydrolase domain-containing protein 9

Immunogen:

Synthesized peptide derived from the Internal region of human ABHD9.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A