

CNTD2 Antibody



PACO05383

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:40000, IHC:1:100-1:300

Protein Background:

Gene ID:

CNTD2

Uniprot

Q9H8S5

Synonyms:

CNTD2; Cyclin N-terminal domain-containing protein 2

Immunogen:

Synthesized peptide derived from the Internal region of human CNTD2.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A