

# CNTNAP3 Antibody



PACO05384

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:10000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

CNTNAP3

**Uniprot**

Q9BZ76

**Synonyms:**

CNTNAP3; CASPR3; KIAA1714; Contactin-associated protein-like 3; Cell recognition molecule Caspr3

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human CASPR3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A