

UBTD1 Antibody



PACO05389

Product Information

Size:

50ug

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

UBTD1

Uniprot

Q9HAC8

Synonyms:

UBTD1; Ubiquitin domain-containing protein 1

Immunogen:

Synthesized peptide derived from the Internal region of human UBTD1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A