

RAB11FIP1 Antibody



PACO05395

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:**Gene ID:**

RAB11FIP1

Uniprot

Q6WKZ4

Synonyms:

RAB11FIP1; RCP; Rab11 family-interacting protein 1; Rab11-FIP1; Rab-coupling protein

Immunogen:

Synthesized peptide derived from the N-terminal region of human Rab11-FIP1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A