

PBX4 Antibody



PACO05400

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:20000, IHC:1:100-1:300

Protein Background:

Gene ID:

PBX4

Uniprot

Q9BYU1

Synonyms:

PBX4; Pre-B-cell leukemia transcription factor 4; Homeobox protein PBX4

Immunogen:

Synthesized peptide derived from the Internal region of human Pbx 4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A