

APOL4 Antibody



PACO05407

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:20000, WB:1:500-1:2000,
IHC:1:100-1:300, IF:1:200-1:1000

Protein Background:**Gene ID:**

APOL4

Uniprot

Q9BPW4

Synonyms:

APOL4; Apolipoprotein L4; Apolipoprotein L-IV; ApoL-IV

Immunogen:

Synthesized peptide derived from the C-terminal region of human ApoL4.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A