

# GPR68 Antibody



PACO05412

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

GPR68

**Uniprot**

Q15743

**Synonyms:**

GPR68; OGR1; Ovarian cancer G-protein coupled receptor 1; OGR-1; G-protein coupled receptor 68; GPR12A; Sphingosylphosphorylcholine receptor

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human GPR68.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A