

# OR51G2 Antibody



PACO05413

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR51G2

**Uniprot**

Q8NGK0

**Synonyms:**

OR51G2; Olfactory receptor 51G2; Olfactory receptor OR11-28

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 51G2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A