

OR2AE1 Antibody



PACO05415

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:10000, IF:1:200-1:1000

Protein Background:

Gene ID:

OR2AE1

Uniprot

Q8NHA4

Synonyms:

OR2AE1; OR2AE2; Olfactory receptor 2AE1; Olfactory receptor 2AE2

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2AE1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A