

# OR2G2 Antibody



PACO05419

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:**

**Gene ID:**

OR2G2

**Uniprot**

Q8NGZ5

**Synonyms:**

OR2G2; Olfactory receptor 2G2; Olfactory receptor OR1-32

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2G2.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A