

OR5V1 Antibody



PACO05426

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

OR5V1

Uniprot

Q9UGF6/Q5SQI9

Synonyms:

Olfactory receptor 5V1

Immunogen:

Synthesized peptide derived from the C-terminal region of human Olfactory receptor 5V1.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A