

# CHST9 Antibody



PACO05462

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300, IF:1:200-1:1000

**Protein Background:****Gene ID:**

CHST9

**Uniprot**

Q7L1S5

**Synonyms:**

CHST9; Carbohydrate sulfotransferase 9; GalNAc-4-O-sulfotransferase 2; GalNAc-4-ST2; GalNAc4ST-2; N-acetylgalactosamine-4-O-sulfotransferase 2

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human CHST9.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A