

DYNLRB2 Antibody



PACO05466

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:10000, IHC:1:100-1:300

Protein Background:**Gene ID:**

DYNLRB2

Uniprot

Q8TF09

Synonyms:

DYNLRB2; DNCL2B; DNLC2B; ROBLD2; Dynein light chain roadblock-type 2; Dynein light chain 2B; cytoplasmic; Roadblock domain-containing protein 2

Immunogen:

Synthesized peptide derived from the Internal region of human Dynein LC 2B.

Storage:

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A