

TAS1R3 Antibody



PACO05475

Product Information

Size: 50ug	Protein Background:
Reactivity: Human	Gene ID: TAS1R3
Source: Rabbit	Uniprot Q7RTX0
Isotype: IgG	Synonyms: TAS1R3; T1R3; TR3; Taste receptor type 1 member 3; Sweet taste receptor T1R3
Applications: ELISA, WB, IF	Immunogen: Synthesized peptide derived from the Internal region of human T1R3.
Recommended dilutions: ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000	Storage: Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Product Images

N/A

N/A