

# OR1D4/OR1D5 Antibody



PACO05480

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:20000, WB:1:500-1:2000

**Protein Background:**

**Gene ID:**

OR1D4/OR1D5

**Uniprot**

P47884/P58170

**Synonyms:**

OR1D4; Olfactory receptor 1D4; Olfactory receptor 17-30; OR17-30; OR1D5; Olfactory receptor 1D5; Olfactory receptor 17-31; OR17-31

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 1D4/1D5.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A