

# OR3A3/OR3A2 Antibody



PACO05484

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:5000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

OR3A3/OR3A2

**Uniprot**

P47888/P47893

**Synonyms:**

OR3A3; OR3A6; OR3A7; OR3A8P; Olfactory receptor 3A3; Olfactory receptor 17-201; OR17-201; Olfactory receptor 3A6; Olfactory receptor 3A7; Olfactory receptor 3A8; Olfactory receptor OR17-22; OR3A2; OLFRA04; Olfactory receptor 3A2; Olfactory

**Immunogen:**

Synthesized peptide derived from the Internal region of human Olfactory receptor 3A2/3A3.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A