

# ULK1 Antibody



PACO05491

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:40000, IHC:1:100-1:300

**Protein Background:**

**Gene ID:**

ULK1

**Uniprot**

O75385

**Synonyms:**

ULK1; KIAA0722; Serine/threonine-protein kinase ULK1; Autophagy-related protein 1 homolog; ATG1; hATG1; Unc-51-like kinase 1

**Immunogen:**

Synthesized peptide derived from the C-terminal region of human ULK1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A