

# QRFPR Antibody



PACO05494

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000, IF:1:200-1:1000

**Protein Background:****Gene ID:**

QRFPR

**Uniprot**

Q96P65

**Synonyms:**

QRFPR; GPR103; Pyroglutamylated RFamide peptide receptor; AQ27; G-protein coupled receptor 103; Orexigenic neuropeptide QRFP receptor; SP9155

**Immunogen:**

Synthesized peptide derived from the Internal region of human GPR103.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A