

# MAF1 Antibody



PACO05507

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:20000, IHC:1:100-1:300

**Protein Background:****Gene ID:**

MAF1

**Uniprot**

Q9H063

**Synonyms:**

MAF1; Repressor of RNA polymerase III transcription MAF1 homolog

**Immunogen:**

Synthesized peptide derived from the Internal region of human Maf1.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A