

# TMEM101 Antibody



PACO05512

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:10000, WB:1:500-1:2000

**Protein Background:****Gene ID:**

TMEM101

**Uniprot**

Q96IK0

**Synonyms:**

TMEM101; Transmembrane protein 101; Putative NF- $\kappa$ B-activating protein 130

**Immunogen:**

Synthesized peptide derived from the Internal region of human TMEM101.

**Storage:**

Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

## Product Images

---

N/A

N/A